

## **INDO-1 LOADING AND SAMPLE STAINING PROCEDURE FOR SIMULTANEOUS MEASUREMENT OF INTRACELLULAR $Ca^{2+}$ AND CELL SURFACE ANTIGEN EXPRESSION**

### **MATERIALS:**

1. 50  $\mu$ g vial Indo-1 (e.g., from LifeTechnologies, Grand Island, NY)
2. DMSO (Sigma-Aldrich, St. Louis, MO)
3. RPMI 1640
4. Monoclonal antibodies (mAb), conjugated to suitable fluorochromes
5. Ionomycin (LifeTechnologies)
6. 37°C water bath, centrifuge, vortexer.
7. Agonists to test  $Ca^{2+}$  flux, e.g. anti-CD3, anti-IgG, ConA.
8. Serum (for RPMI with 2% serum, if cells require serum).

### **METHOD:**

1. Incubate cells ( $\leq 2 \times 10^7$ /mL) in RPMI with 1-5  $\mu$ M Indo-1 (acetoxymethyl ester) at 37°C for 40 min for loading.
2. Incubate aliquots of Indo-1 loaded cells with saturating concentrations of fluorochrome-conjugated antibodies for 20 min. Incubate at 20° to 25°C unless the antigen is subject to capping, otherwise use 4° to 8°C. Note: mAbs must be azide free. Note: set-up single color stained cells for setting appropriate fluorescence compensation on the instrument.
3. Wash cells twice in RPMI and suspend them at the desired concentration (usually  $2 \times 10^6$ /mL). Higher cell concentrations ( $4 \times 10^6$ /mL) are required when the cells of interest represent less than 10% of the total population. Cells can be kept at 20° to 25°C unless the antigen is subject to capping, otherwise use 4° to 8°C.
4. Samples should be analyzed shortly after the cells were prepared.
5. Ionomycin (1-3  $\mu$ M final conc.) is used as a positive control for Indo-1 loading and maximum  $Ca^{2+}$  flux.

## PREPARATION OF INDO-1:

1. Add 150  $\mu\text{L}$  of DMSO to a 50  $\mu\text{g}$  vial of Indo-1, cover with aluminum foil to protect from light.
2. Vortex well, then warm to 37°C for 5 min.
3. Transfer 150  $\mu\text{L}$  of Indo-1 from vial to 4.85 mls of RPMI (=10 $\mu\text{M}$ ). Wash out vial very well. If not the entire amount of Indo-1 dissolved in DMSO is used, the remainder can be stored dessicated at -20°C for a maximum of 6 months.
4. Cover the tube of 5 mL of 10  $\mu\text{M}$  Indo-1 with foil.
5. Aliquot the appropriate amount of 10 $\mu\text{M}$  Indo-1 to the cell suspension (final conc.=1-5  $\mu\text{M}$ ). The optimal concentration is dependent on the cell type.
6. Store excess RPMI-diluted 10  $\mu\text{M}$  Indo-1 at 4°C. In our laboratory, the 10  $\mu\text{M}$  Indo-1 solution has been tested for stability up to 24 hrs.

## PREPARATION OF IONOMYCIN:

1. Dissolve 1 mg of ionomycin in 1 mL DMSO.
2. Aliquot 13.5  $\mu\text{L}$  of ionomycin solution into vials for later use and store at -20°C for less than one year.
3. Dilute one 13.5 $\mu\text{L}$  vial of ionomycin with RPMI to a volume of 3 mL (=6  $\mu\text{M}$ ).
4. Cover the 3 mL of 6  $\mu\text{M}$  working stock with foil to protect from light.
5. 150  $\mu\text{L}$  of working stock ionomycin is added to 300  $\mu\text{L}$  of Indo-1 loaded cell suspension.

### Special Note:

Indo-1 requires UV excitation and special filter sets for optimal detection, e.g., a 395-415 nm bandpass (BP) filter for Indo 'violet' and a 515-545 nm BP for Indo 'blue-green'.

### Further Reading:

June CH, Rabinovitch PS. Intracellular ionized calcium. *Methods Cell Biol.* 1994;41:149-74.

June CH, Abe R, Rabinovitch PS. Measurement of intracellular calcium ions by flow cytometry. In: Current Protocols in Cytometry, Vol 2, Robinson JP, Darzynkiewicz Z, Hyun W, Orfao A, Rabinovitch P, eds., John Wiley & Sons, 1997, pp. 9.8.1 – 9.8.19.