

BLOCKING OF FC-RECEPTOR BINDING WITH HUMAN AB SERUM (HAB)

If your cells have many Fc-receptors on the cell surface (in particular monocytes, macrophages) or they have been cultured in serum-free medium, it is advisable to block nonspecific binding of monoclonal antibodies by pre-incubation of cells with human AB serum (HAB). Note that for staining of whole blood this is not necessary, because serum in high concentration is present during staining.

I. REAGENTS:

A. Antibodies

1. Primary antibodies: usually purchased or your own monoclonals

a. If these are conjugated with a fluorochrome (i.e., FITC, PE, or other) use the DIRECT STAINING procedure.

b) If the primary antibody is not conjugated to a fluorochrome, use the INDIRECT STAINING procedure.

2. Secondary antibodies: fluorochrome-conjugated polyclonal antibody

B. Buffer: PBS (Ca^{2+} and Mg^{2+} free, e.g., Cat. #9240, Irvine Scientific, CA) + 2% newborn calf serum (or 0.2% BSA) + 0.1% sodium azide.

C. HAB (e.g., from Innovative Research, MI), either purchase heat-inactivated HAB or inactivate it by heating it to 56°C for 1h. Divide into small aliquots and keep frozen at -20°C .

II. PROCEDURE:

DIRECT STAINING

1. Prepare your cells as a suspension of single cells in a manner appropriate for the specimen you wish to examine. As the final step, wash at least once with 1 ml of cold BUFFER. Re-suspend the cells at 10^7 cells/ml (thus $50\ \mu\text{L} = 5 \times 10^5$ cells) in cold BUFFER.

Check cell viability, it should exceed 90%. If cell viability is less than 90%, remove dead cells by Ficoll-Hypaque separation, otherwise dead cells will bind antibodies nonspecifically. Preferably, add a dead cell discriminating dye for exclusion of dead cells from flow cytometric analysis (see protocols for adding either propidium iodide or 7-amino-actinomycin D to the cells in the final re-suspension step before acquisition on the flow cytometer or follow procedures as outlined by manufacturers of other dead cell discrimination dyes).

2. Add $50\ \mu\text{L}$ of the cell suspension to the bottom of the tubes.

3. Add $50\ \mu\text{L}$ of HAB to each tube mix well and incubate for ~ 1 min at room temperature.

4. Then, add the appropriate amount of monoclonal antibody to the bottom of tubes sitting on ice.

Note: For multi-color staining, add all your fluorochrome-conjugated antibodies at the same time.

5. Vortex briefly and incubate for 30 min at 4°C in the dark.
6. Wash twice with 1mL of buffer and hold them at 4°C (or on ice) prior to acquisition on the flow cytometer.

INDIRECT STAINING

- 1-6. Process cell samples as above using a working dilution of unlabelled monoclonal antibody.
7. Resuspend cell pellet in 50 µL HAB, mix and incubate for ~1 min at room temperature.
8. Add 50 µL of working dilution of the fluorochrome-conjugated second antibody.
9. Vortex briefly and incubate for 20 min at 4°C in the dark.
10. Wash twice with 1 mL of buffer; centrifuge at 250g for 5 min.
11. Resuspend samples in 1 mL of buffer and hold them at 4°C (or on ice) protected from light prior to acquisition on the flow cytometer.

Note: This method cannot be used for staining of surface Ig or staining with antibodies that are directed against Fc-receptors (e.g., CD16).